


COUNT IT DOWN

Foley Catheter Deflation


1 Exercise the syringe

Slide the syringe plunger up and down the barrel several times to "loosen it up."



2 Back-off to 0.5cc


Compress the plunger all the way and then pull back 0.5cc to prevent adherence.



3 Let balloon volunteer its water

Gently insert syringe in catheter inflation valve, using only enough force to make the syringe "stick" in the valve.

Allow pressure within the balloon to force the plunger back as it fills with water.




Do not pull back on plunger!

4 Wait at least 30 seconds

It takes thirty seconds for the water to flow through the small-diameter inflation lumen.

If you notice slow or no deflation, gently re-seat the syringe.



IF NO DEFLATION



5 Reposition patient


If the retention balloon still does not deflate, reposition the patient to ensure catheter is not in traction or compressed within the bladder, and urine flows freely.



If balloon still fails to deflate, apply gentle, slow aspiration. Forceful or rapid aspiration may cause the inflation lumen within the Foley catheter to collapse.

If problems still persist after all previous steps have been followed, refer to hospital protocol and contact an adequately trained professional for assistance as directed.


For clinical consultation, available 24-hours a day, 7 days a week, contact the BARD® Medical Services and Support Team at 800.227.3357.



5cc Foley catheter inflated with 10cc of water

FILL IT UP

INFLATION REMINDERS



5cc Foley catheter inflated with 5cc of water

Assure catheter balloon is positioned well-within the patient's bladder.

Balloon should be inflated slowly with a gentle, constant force.

Completely inflate 5cc balloon with the volume prescribed on the package using the entire 10cc of sterile water provided.

Improperly inflated Foley catheter balloons may cause drainage and deflation difficulties

BARD | MEDICAL

Advancing Lives and the Delivery of Healthcare.™

C. R. Bard, Inc.
8195 Industrial Blvd
Covington, GA 30014
1-800-526-4455
www.bardmedical.com

Foley catheters are intended for use in the drainage and/or collection and/or measurement of urine.

Cautions:

- Latex Foley catheters contain natural rubber latex which may cause allergic reactions.
- Federal (U.S.A.) law restricts these devices to sale by or on the order of a physician.

Warnings:

- On catheter, do not use ointments or lubricants having a petrolatum base. They will damage the catheter and may cause balloon to burst.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practices and applicable local, state and federal laws and regulations.

Please consult product label and insert for any indications, contraindications, hazards, warnings, cautions and directions for use.

Bard and Advancing Lives and the Delivery of Healthcare are trademarks and/or registered trademark of C. R. Bard, Inc. ©2017 C. R. Bard, Inc. All rights reserved. 1706-32 BMD/FOCA/0717/0053