




For Capture of Challenging Calculi...



Bard: Ureteral Stent

---

**Bard® Dimension®**  
ARTICULATING STONE BASKET/GRASPER


The complex block features a decorative header with a blue and green geometric pattern. Below this, the text "Bard: Ureteral Stent" is centered. A horizontal line separates this from the main product name, "Bard® Dimension®", which is rendered in a large, bold, teal font. Underneath the product name, the text "ARTICULATING STONE BASKET/GRASPER" is centered in a smaller, teal, all-caps font.



## For Capture of Challenging Calculi...

### Articulating feature combines the advantages of a basket and grasper in a single instrument

To help you meet the challenge presented by calculi in difficult locations, we offer an innovative new solution: the DIMENSION® Articulating Stone Basket/Grasper. An impressive combination of design features allows you to capture and release calculi and fragments with enhanced efficiency and minimal risk of trauma. We are confident this revolutionary technology will add a valuable new dimension of control in a variety of stone removal procedures.




\*Data on file, C. R. Bard, Inc.



## For Capture of Challenging Calculi...

### Efficient capture and controlled release

- New articulating feature offers the advantages of a basket and grasper in a single instrument.
- Unique “butterfly” wire configuration creates a tipless basket - for easier capture, controlled release and minimal risk of trauma.
- Teardrop shape permits efficient capture of a wide range of stones; Nitinol wire is kink-resistant and retains its shape even in tortuous anatomy.



### Engineered for precision control

- Bard® Platinum Class™ handle helps assure ergonomic comfort and single-hand control of either passive or articulating capability.
- Torque-stable sheath and drive wire promote effective navigation through tortuous anatomy.






## For Capture of Challenging Calculi...

### Achieves closer proximity to stone

- Proprietary “butterfly” design permits closer proximity to calculi than other widely used baskets.
- Unlike conventional designs, offers ability to grasp and engage calculi directly in front of endoscope, allowing easier capture of stones even in lower pole of kidney.



### Delivers superior performance

- In a laboratory test using a simulated stone,\* the DIMENSION® basket demonstrated an 83% better capture rate than a leading competitor on the first capture attempt.
- In up to 3 attempts, the DIMENSION® basket demonstrated a 73% better capture rate than a leading competitor.



\*Data on file, C. R. Bard, Inc.



For Capture of Challenging Calculi...

## Precision Control. Proven Performance.

Order #	Balloon Diameter (mm)	Length (cm)	Shaft Size (Fr)
042310	10	115	2.4
042313	13	115	2.4
042316	16	115	2.4
043310	10	115	3
043313	13	115	3
043316	16	115	3



Where Knowledge Leads™

C. R. Bard, Inc. Covington, GA 30014 888.367.2273  
www.bardurological.com



Bard and Dimension are registered trademarks of C. R. Bard, Inc. or an affiliate. Where Knowledge Leads and Take a Closer Look are trademarks of C. R. Bard, Inc. or an affiliate.

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use. © 2006 C. R. Bard, Inc. All Rights Reserved. Printed in U.S.A. 0708-37 R8/07