



Products Designed for Procedural Flexibility



Ureteral Catheters

Designed for Procedural Flexibility


Bard offers a full range of ureteral catheters designed to meet any procedural need. With a variety of materials, tip configurations and sizes, the Bard Ureteral Catheter choices are specially engineered for ease of navigation with minimal risk of tissue trauma. Convenient shaft color coding aids in French size identification.

TigerTail® Ureteral Catheter



- Soft, smooth and flexible tip tracks easily and efficiently through difficult anatomy.

Polyurethane Ureteral Catheters (disposable)




- High quality, single use polyurethane ureteral catheters that provide excellent handling characteristics, flexibility and resistance to kinking.

Woven Ureteral Catheters (re-sterilizable)



- Premium quality catheters that feature a tightly woven Dacron® core to prevent kinking and provide a large, smooth drainage lumen

Dual Lumen Ureteral Catheter



- Designed with two independent lumens for over-the-wire placement and simultaneous fluid injection or secondary wire placement.

Ureteral Catheters

Products Designed for Procedural Flexibility

Ordering Information

Woven Ureteral Catheter

Product Description	Catalog#	Size	Number of Units
Whistle Tip	036103	70cm, 3 FR	1
	036104	70cm, 4 FR	1
	036105	70cm, 5 FR	1
	036106	70cm, 6 FR	1
Round Tip	036303	70cm, 3 FR	1
	036304	70cm, 4 FR	1
	036305	70cm, 5 FR	1
	036306	70cm, 6 FR	1
Olive Tip	036503	70cm, 3 FR	1
	036504	70cm, 4 FR	1
	036505	70cm, 5 FR	1
	036506	70cm, 6 FR	1
Blasucci Flexible Filiform Tip	037404	70cm, 4 FR	1
	037405	70cm, 5 FR	1
	037406	70cm, 6 FR	1
Blasucci Spiral Filiform Tip	037504	70cm, 4 FR	1
	037505	70cm, 5 FR	1
	037506	70cm, 6 FR	1
	037508	70cm, 8 FR	1
Braasch Bulb Whistle Tip	037606	70cm, 4 FR	1
	037607	70cm, 5 FR	1
	037608	70cm, 6 FR	1
	037610	70cm, 10 FR	1
Cone Tip	038008	70cm, 8 FR	1
	038010	70cm, 10 FR	1
	038012	70cm, 12 FR	1
Wishard Flat Coude Tip	038604	70cm, 4 FR	1
	038605	70cm, 5 FR	1
	038606	70cm, 6 FR	1
Braasch Bulb Terminal Eye Only	603708	70cm, 8 FR	1
	603710	70cm, 10 FR	1

Polyurethane Ureteral Catheter

Product Description	Catalog#	Guidewire Compatibility	Size	Number of Units
Open Tip without side eyes	134005	.038"	69cm, 5 FR	10
	134006	.038"	69cm, 6 FR	10
	134104	.025"	70cm, 4 FR	10
	134105	.038"	70cm, 5 FR	10
	134106	.038"	70cm, 6 FR	10
	with side eyes	134206	.038"	70cm, 6 FR
135005		.038"	69cm, 5 FR	10
135006		.038"	69cm, 6 FR	10
135105		.038"	115cm, 5 FR	10
135106		.038"	115cm, 6 FR	10
Whistle Tip		136404		70cm, 4 FR
	136405		70cm, 5 FR	10
	136406		70cm, 6 FR	10
Olive Tip	136504		70cm, 4 FR	10
	136505		70cm, 5 FR	10
	136506		70cm, 6 FR	10
Spiral Tip	137504		70cm, 4 FR	10
	137505		70cm, 5 FR	10
	137506		70cm, 6 FR	10
Cone Tip	138008		70cm, 5 FR	10
	138010		70cm, 6 FR	10
	138012		70cm, 7 FR	10

TigerTail® Ureteral Catheter

Product Description	Catalog#	Guidewire Compatibility	Size	Number of Units
Open Tip, without side eyes	139004	.025"	70cm, 4 FR	10
	139005	.038"	70cm, 5 FR	10
	139006	.038"	70cm, 6 FR	10
Open Tip, with Side Eyes	139104	.025"	70cm, 4 FR	10
	139105	.038"	70cm, 5 FR	10
	139106	.038"	70cm, 6 FR	10

Product Description	Catalog#
Dispenser Wall Rack	1369
Ureteral Catheter Plastic Adapter	140000
Heyman™ Slipover Catheter Stylet, 6 Fr.	042900

Product Description	Catalog#	Guidewire Compatibility	Size	Number of Units
Dual Lumen Ureteral Catheter	130200	.038"	50cm, 10 FR	1

For additional information contact your local Bard representative or call 1.888.367.2273.



Where Knowledge Leads™

C. R. Bard, Inc. Covington, GA 30014

www.bardurological.com



Bard and TigerTail are registered trademarks of C. R. Bard, Inc. or an affiliate. Where Knowledge Leads is a trademark of C. R. Bard, Inc. or an affiliate. Dacron is a registered trademark of Invista.
© 2006 C. R. Bard, Inc. All Rights Reserved. Printed in U.S.A.

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

0606-15 R06/06 THP P07/06 5M