

INTERMITTENT CATHETERS

CLEAN-CATH[®]

INTERMITTENT CATHETER

The Clinicians Choice for
Self Catheterization



UROLOGY

BARD

CLEAN-CATH[®]

INTERMITTENT CATHETER

The Clinicians Choice
for Self Catheterization

Intermittent Catheter

- Finger grips for ease of use
- Flexible funnel
- Available with convenient color-coded funnels
- Latex Free
- Individually packaged in easy to open peel-strip package
- Sterile, Single-Use

Straight Tip

- Range of French Sizes, 6 Fr. - 18 Fr.
- Smooth, round tip for easy insertion

Coude Tip

- Range of French Sizes, 10 Fr. - 18 Fr.
- Olive Coude tip catheter makes it easier to thread the catheter past the prostatic curve
- Coude tip catheters are also helpful for insertion difficulties due to narrowing of the urethra and blockages in the urethra pathway

Regular Tip Ordering Information

Reorder Code	French Size	Length (inches)	Funnel Color	Case Quantity
420706	6	6	Green	50
420708	8	6	Blue	50
420710	10	6	Black	50
420712	12	6	White	50
420714	14	6	Green	50
420716	16	6	Orange	50
420718	18	6	Red	50
421108	8	10	Blue	50
421110	10	10	Black	50
421708	8	16	Blue	50
421710	10	16	Black	50
421712	12	16	White	50
421714	14	16	Green	50
421716	16	16	Orange	50
421718	18	16	Red	50

Coude Tip Ordering Information

Reorder Code	French Size	Length (inches)	Funnel Color	Case Quantity
423710	10	16	Black	30
423712	12	16	White	30
423714	14	16	Green	30
423716	16	16	Orange	30
423718	18	16	Red	30

For more information about the CLEAN-CATH[®] Intermittent Catheter and other Bard products, please contact your Bard sales representative, call 800.526.4455, or visit www.bardmedical.com.

BARD

C. R. Bard, Inc.
Covington, GA 30014-2655

Bard and Clean-Cath are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and directions for use.

©2012 C. R. Bard, Inc. All Rights Reserved. 1201-69 R03/12 THP P03/12 5M

BARD

INTERMITTENT CATHETERS

INTERMITTENT CATHETERS