

Bard® Foley Catheter Selection Guide

Catalog No.*	Category	Balloon Volume	Available Sizes (Fr.)**	Material	Coating	Tip Length
0100LXX	2-Way	30cc	16-24	Amber Latex	Bard® Hydrogel	Med.
0102LXX	2-Way	5cc	12-24	Red Latex	Bard® Hydrogel	Med.
0102SIXX	2-Way	5cc	12-24	Red Latex	Baci-Guard™ silver alloy with Bard® Hydrogel	Med.
0103LXX	2-Way	30cc	14-24	Red Latex	Bard® Hydrogel	Med.
0103SIXX	2-Way	30cc	14-24	Red Latex	Baci-Guard™ silver alloy with Bard® Hydrogel	Med.
0113LXX	2-Way	75cc	20-28	Amber Latex	Bard® Hydrogel	Med.
0118LXX	2-Way	30cc	16-24	Amber Latex	Bard® Hydrogel	Short
0119LXX	3-Way	5cc	16-26	Amber Latex	Bard® Hydrogel	Short
0119SIXX	3-Way	5cc	16-26	Amber Latex	Baci-Guard™ silver alloy with Bard® Hydrogel	Short
012375LXX	2-Way	75cc	22-26	Amber Latex	Bard® Hydrogel	Med.
0123RLXX	2-Way	30cc	20-26	Red Latex	Bard® Hydrogel	Med.
0125RLXX	3-Way	30cc	20-24	Red Latex	Bard® Hydrogel	Med.
012675XX	3-Way	75cc	20-24	Amber Latex	Bard® Hydrogel	Med.
0132LXX	3-Way	5cc	18-24	Amber Latex	Bard® Hydrogel	Med.
0134LXX	3-Way	30cc	18-26	Amber Latex	Bard® Hydrogel	Med.
0136LXX	2-Way	30cc	22 & 25	Amber Latex	Bard® Hydrogel	Long
0138LXX	3-Way	30cc	18-26	Amber Latex	Bard® Hydrogel	Long
0149RL24	2-Way	30cc	24	Red Latex	Bard® Hydrogel	Short
0167LXX	3-Way	30cc	16-26	Amber Latex	Bard® Hydrogel	Med.
0167SIXX	3-Way	30cc	16-26	Latex	Baci-Guard™ silver alloy with Bard® Hydrogel	Med.
0167VXX	3-Way	30cc	16-26	Amber Latex	Silicone-Elastomer	Med.
0168LXX	2-Way	5cc	12-24	Amber Latex	Bard® Hydrogel	Med.
0168SI	2-Way	5cc	12-24	Amber Latex	Baci-Guard™ silver alloy with Bard® Hydrogel	Med.
0169LXX	2-Way	5cc	12-24	Amber Latex	Bard® Hydrogel	Short
0171L	Diagnostic	5cc	18	Amber Latex	Bard® Hydrogel	Short
0174LXX	3-Way	5cc	18-26	Amber Latex	Bard® Hydrogel	Med.
0187L	Diagnostic	20cc/30cc	14	Amber Latex	Bard® Hydrogel	Short
0193L24	3-Way	30cc	24	Amber Latex	Bard® Hydrogel	X-Long
0196LXX	2-Way	5cc	16-22	Red Latex	Bard® Hydrogel	Short
0196SIXX	2-Way	5cc	16-22	Red Latex	Baci-Guard™ silver alloy with Bard® Hydrogel	Short
0198LXX	2-Way	30cc	22-26	Red Latex	Bard® Hydrogel	Short
1851SIXX	3-Way	30cc	18-24	Red Latex	Baci-Guard™ silver alloy with Bard® Hydrogel	Long
1853SIXX	3-Way	30cc	18-24	Red Latex	Baci-Guard™ silver alloy with Bard® Hydrogel	Med.

Bard® Foley Catheter Selection Guide

Tip Design	Drainage Eyes	Model Name/ Other Features	Page #
Coude	Single Eye	Carson model	20
Coude	Single Eye	Tiemann model	21
Coude	Single Eye	Tiemann model	8
Coude	Single Eye	Tiemann model	21
Coude	Single Eye	Tiemann model	9
Round	2 Slaggered	Extra-strong ovoid balloon	22
Round	2 Opposed		23
Round	2 Opposed	3-Way for continuous irrigation	35
Round	2 Opposed	3-Way for continuous irrigation\	12
Round	2 Slaggered		24
Round	2 Slaggered	Radiopaque	25
Round	2 Slaggered	3-Way for continuous irrigation, Radiopaque	36
Round	2 Slaggered	Creivy Model, long 75cc balloon	37
Open Whistle	2 on same side	3-Way for continuous irrigation	38
Open Whistle	2 on same side	Alcock Model	39
Round	3 Slaggered		26
Round	3 Slaggered	Emmett Model	40
Open Whistle	Single Eye	Ainsworth model, Radiopaque	27
Round	2 Slaggered	3-Way for continuous irrigation	41
Round	2 Slaggered	3-Way for continuous irrigation	13
Round	2 Slaggered	3-Way for continuous irrigation	41
Coude	Single Eye	Carson model	28
Coude	Single Eye	Carson model	10
Round	2 Opposed	Female Length Catheter Shaft (8"	29
Round	2 Opposed	Lapides model,	51
Round	2 Slaggered	3-Way for continuous irrigation	42
Round	One, between balloons	Davis Model	52
Natural Curve	3 Slaggered	Suggs Model	43
Open	2 Opposed	Councill model, Radiopaque	31
Councill	2 Opposed	Councill model, Radiopaque	11
Round	Single, Large Eye	Wheeler model, Radiopaque	32
Round	2+2 Opposed	3-Way for continuous irrigation, Stiff Hemostasis	14
Open Concave	2 Opposed	3-Way for continuous irrigation, Stiff Hemostasis	15

Bard® Foley Catheter Selection Guide

Catalog No.*	Category	Balloon Volume	Available Sizes (Fr.)**	Material	Coating	Tip Length
1857SI	3-Way	30cc	18-24	Red Latex	Back-Guard™ silver alloy with Bard® Hydrogel	Med.
2550HXX	2-Way	30cc	20-24	Amber Latex	Bard® Hydrogel	Med.
2556HXX	2-Way	30cc	20-24	Amber Latex	Bard® Hydrogel	Long
2551HXX	3-Way	30cc	20-24	Amber Latex	Bard® Hydrogel	Long
2557HXX	3-Way	30cc	20-24	Amber Latex	Bard® Hydrogel	Long
334XX	2-Way	30cc	16-28	Latex	Silicone-Elastomer	Med.
336XX	2-Way	5cc	14-24	Latex	Silicone-Elastomer	Short
6003LXX	3-Way	30cc	22 & 24	Red Latex	Bard® Hydrogel	Long
6016L22	Hemostatic	40cc/20cc	22	Red Latex	Bard® Hydrogel	Med.
6016L24	Hemostatic	5cc/15cc	24	Amber Latex	Bard® Hydrogel	Short
600190L	Diagnostic	20cc/5cc	20	Amber Latex	Bard® Hydrogel	Short
6097L24	3-Way	30cc	24	Amber Latex	Bard® Hydrogel	Short
705XXL	3-Way	5cc	16-24	Silicone	Bard® Hydrogel	Med.
705XXSI	3-Way	5cc	16-24	Silicone	Back-Guard™ silver alloy with Bard® Hydrogel	Med.
730XXL	3-Way	30cc	16-24	Silicone	Bard® Hydrogel	Med.
730XXSI	3-Way	30cc	16-24	Silicone	Back-Guard™ silver alloy with Bard® Hydrogel	Med.

*Substitute the balloon size for the "xx" in the catalog number designation
 (for 8 Fr. catheter shafts, substitute "08") **Even French (Fr.) sizes unless otherwise noted

Bard® Foley Catheter Selection Guide

Tip Design	Drainage Eyes	Model Name/ Other Features	Page #
Open Concave	2-Staggered	Couvelaire	16
Round	2 sets of 2 Opposed	Hematuria, Coil reinforced shaft	33
Open Coude	2 Staggered	Hematuria, Coil reinforced shaft	34
Round	2 sets of 2 Opposed	Hematuria, Coil reinforced shaft	44
Open Coude	2 Staggered	Hematuria, Coil reinforced shaft	45
Round	2 Staggered	SILASTIC®	30
Round	2 Opposed	SILASTIC®	30
Coude	Single Eye	Hotter Model, Radiopaque	46
Open Whistle	2 on same side	Coleman model	47
Open Whistle	Single Eye	Lerman model	48
Round	One, between balloons	Trattner Model	53
Round	2 Opposed	Naegeli Model	49
Round	2-Opposed	3-Way, Latex-Free, Lubri-Sil®	50
Round	2-Opposed	3-Way, Latex-Free, Lubri-Sil® IC	17
Round	2-Opposed	3-Way, Latex-Free, Lubri-Sil®	50
Round	2-Opposed	3-Way, Latex-Free, Lubri-Sil® IC	17