


TIGERTAIL®

URETERAL CATHETER



Flexible Open Tip Polyurethane Ureteral Catheter

- Ideal for navigating tortuous pathways
- Soft, smooth and flexible tip tracks easily and efficiently through tortuous anatomy
- Helps promote atraumatic ureteral access
- Radiopaque markers at distal end permit clear visualization under fluoroscopy


Warning: Do not use as a ureteral splint.

Please consult product inserts and labels for any indications, contraindications, hazards, warnings, cautions and directions for use.

Bard and TigerTail are trademarks and/or registered trademarks of C. R. Bard, Inc.

©2011 C. R. Bard, Inc. All Rights Reserved. Printed in the U.S.A.

1105-06 R05/11 P05/11



Cook Pollack Open-End Flexi-Tip® Ureteral Catheter BARD® TIGERTAIL® Flexible Open Tip Ureteral Catheter Microvasive Flexible Tip Open-End Ureteral Catheter

Fluoroscopic Comparison of BARD® TIGERTAIL® Versus Catheters from Two Competitors.

Product Description	Catalog #	Guidewire Compatibility	Size	Number of Units
Open Tip without side eyes	139004	.025"	70cm, 4 Fr	10
	139005	.038"	70cm, 5 Fr	10
	139006	.038"	70cm, 6 Fr	10
Open Tip, with side eyes	139104	.025"	70cm, 4 Fr	10
	139105	.038"	70cm, 5 Fr	10
	139106	.038"	70cm, 6 Fr	10