

Bard® Clean-Cath® Intermittent Catheter

Bard Clean-Cath Intermittent Catheter

The Clinicians Choice for Self Catheterization

Complete Range of French Sizes

6 Fr. - 18 Fr.

Color Coded

Available with convenient color-coded funnels

Easy to Use

Smooth, round tip for easy insertion

Latex Free

100% Latex Free


Individually Packaged

Individually packaged in easy to open peel-strip package

Sterile, Single-Use

Phone Number:

1-800-526-4455 or contact your Bard Representative



Ordering Information

<i>Reorder Code</i>	<i>French Size</i>	<i>Length (Inches)</i>	<i>Funnel Color</i>	<i>Case Quantity</i>
420706	6	6	Green	50
420708	8	6	Blue	50
420710	10	6	Black	50
420712	12	6	White	50
420714	14	6	Green	50
420716	16	6	Orange	50
420718	18	6	Red	50
421108	8	10	Blue	50
421110	10	10	Black	50
421708	8	16	Blue	50
421710	10	16	Black	50
421712	12	16	White	50
421714	14	16	Green	50
421716	16	16	Orange	50
421718	18	16	Red	50

Bard Medical Division

C. R. Bard, Inc.

Covington, GA 30014-2655

Bard and Clean-Cath are registered trademarks of C. R. Bard, Inc. or an affiliate. Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and directions for use. © 2002 C. R. Bard, Inc. All rights reserved. 0112-08 P01/02 THP P01/02 54